

## The Hyper-Realistic Moulage Kit (HRMK)

The Hyper-Realistic Moulage Kit provides the capabilities of a professional Special Effects Make-up Artist in one highly versatile product. The HRMK inserts Hyper-Realism into any training revolution, increasing the stress inoculation factor for medical and non-medical professionals; greatly decreases the costs when implementing Hyper-Realistic™ casualties into training scenarios because everything is completely reusable. The HRMK is 100% user customizable enabling the order of specific "pages" with the desired wound patterns or types.



Custom Designed and Fabricated  
Carrying Case



### The HRMK contains:

- All required supplies needed for wound application and removal
- The ability to produce up to 12.5 gallons of Artificial Blood
- Step by step instructions for application, removal, cleaning, reusing, and storing of wounds
- Reusable and customizable injury reference cards to increase the realism for casualties as well as simplify wound description and application for non medical personnel
- Removable "pages" of wound patterns for ease of use and accountability

### Customizable wound patterns:



- Superficial and deep lacerations
- Closed and compound fractures
- Impaled objects and shrapnel wounds
- Partial and full thickness burns
- Partial upper amputations
- Abdominal eviscerations
- Flail chest and open thoracic injuries
- Extraorbital avulsions with or without ocular injury
- Maxiofacial trauma with and without dental injuries
- Chemical, Biological, Radiological and Nuclear injuries
- Entrance and exit gunshot wounds from multiple calibers
  - 5.56, 7.62, and 9mm

**All wounds can be applied on any anatomical location to create the desired effect**

### Unique Features:

- Complete re-usability of wounds
- All supplies stored in one portable bag for ease of use and convenience

**The Hyper-Realistic Moulage Kit (HRMK)**

|  |   |  |
|--|---|--|
| <p>A1</p> <p><b>STRATEGIC OPERATIONS</b><br/> <small>HYPER-REALISTIC TACTICAL TRAINING THE NEXT GENERATION</small></p> <p><u>Wound Description</u><br/>                 Compound Ulnar fracture<br/>                 Deep forearm laceration</p> <p><u>Signs &amp; Symptoms</u><br/>                 Severe pain<br/>                 Signs of shock</p> |  |  |
|--|---|--|

(Front)

Injury Reference Card for each Wound

(Back)



Eye Avulsion



Compound Humerous Fracture



Dorsal Hand 1st & 2nd Degree Burns



Superficial Large Laceration



Bowel Evisceration



Gunshot Entrance Wound



Forearm Laceration with Shrapnel



Upper Extremity 3rd Degree Burn



Compound Fibula Fracture



Gunshot Exit Wound



Penetrating Shrapnel Wounds



Sucking Chest Wound



2nd Degree Partial Thickness Burn



Forearm Laceration



Partial Hand Amputation



Partial Foot Amputation